



Price list

valid from 01.01.2024

| Product ID | Description | Price / unit |
|--------------------|--|--------------|
| | | |
| 3D Skin Models | | |
| FT SM-1 | Phenion® Full-Thickness Skin Model (single tissue), ø 1.4 cm | \$ 168 |
| FT INSERT-1 | Phenion® Full-Thickness INSERT Skin Model (single tissue), ø 1.4 cm | \$ 168 |
| FT LARGE-1 | Phenion® Full-Thickness LARGE Skin Model (single tissue), ø 3.1 cm | \$ 282 |
| FT SMALL-1 | Phenion® Full-Thickness SMALL Skin Model (single tissue), Ø 0.9 cm, on filter or insert-based | \$ 167 |
| FT LL-1 | Phenion® Full-Thickness LONG-LIFE Skin Model (single tissue), ø 1.4 cm | \$ 228 |
| FT AGED-1 | Phenion® Full-Thickness AGED Skin Model (single tissue), ø 1.4 cm | \$ 180 |
| Human Skin Cells | | |
| hFB P1 | Phenion® Fibroblasts, human foreskin, juvenile P1 (1x10*6) | \$ 600 |
| hK P1 | Phenion® Keratinocytes, human foreskin, juvenile P1 (1x10*6) | \$ 697 |
| hFeeder | Phenion® Feeder Cells, derived from human foreskin fibroblasts (1x10*6) | \$ 697 |
| Cell Culture Media | | |
| FB CM-250 | Phenion® Fibroblast Culture Medium (250 mL) | \$ 66 |
| K CM-250 | Phenion® Keratinocyte Culture Medium (250 mL) | \$ 102 |
| ALI CM-250 | Phenion® Air-Liquid Interface Culture Medium (250 mL) | \$ 90 |
| ALI CM PR-250 | Phenion® Air-Liquid Interface Culture Medium w/o phenol-red (250 mL) | \$ 132 |
| ALI CM HC-250 | Phenion® Air-Liquid Interface Culture Medium w/o hydrocortisone (250 mL) | \$ 132 |
| ALI CM AB-250 | Phenion® Air-Liquid Interface Culture Medium w/o antibiotics (250 mL) | \$ 132 |
| Accessories | | |
| FCA-15 | Phenion® Franz Cell Adapter The tailored solution for skin permeation experiments with Phenion® FT Skin Models (fits to 15 mm PermeGear Franz Cells) | \$ 180 |

We allow discounts in case of exceeding following order volumes: > \$ 3.600 = 5%, > \$ 7.200 = 7%, > \$ 12.000 = 9%





