

## Phenion product portfolio

Product-ID	Description	Price / unit
<b>3D Skin Models</b>		
FT SM-1	Phenion® Full-Thickness Skin Model (single tissue), ø 1.4 cm	\$ 150
FT INSERT-1	Phenion® Full-Thickness INSERT Skin Model (single tissue), ø 1.4 cm	\$ 150
FT LARGE-1	Phenion® Full-Thickness LARGE Skin Model (single tissue), ø 3.1 cm	\$ 282
FT SMALL-1	Phenion® Full-Thickness SMALL Skin Model (single tissue), ø 0.9 cm, also available in insert	\$ 162
FT LL-1	Phenion® Full-Thickness LONG-LIFE Skin Model (single tissue), ø 1.4 cm	\$ 228
FT AGED-1	Phenion® Full-Thickness AGED Skin Model (single tissue), ø 1.4 cm	\$ 180
CS-1001	<b>NEW epiCS®</b> Epidermal Model (single tissue), ø 0.9 cm	\$ 63
CS-1101	<b>NEW epiCS®-M</b> Epidermal Model w/ melanocytes, Caucasian donor, (single tissue), ø 0.9 cm	\$ 108
CS-1111	<b>NEW epiCS®-M</b> Epidermal Model w/ melanocytes, Asian-Caucasian donor, (single tissue), ø 0.9 cm	\$ 108
CS-1121	<b>NEW epiCS®-M</b> Epidermal Model w/ melanocytes, Afro-American donor, (single tissue), ø 0.9 cm	\$ 108
<b>Human Skin Cells</b>		
hFB P1	Phenion® Fibroblasts, human foreskin, juvenile P1 (1x10*6)	\$ 600
hK P1	Phenion® Keratinocytes, human foreskin, juvenile P1 (1x10*6)	\$ 697
hFeeder	Phenion® Feeder Cells, derived from human foreskin fibroblasts (1x10*6)	\$ 697
<b>Cell Culture Media</b>		
FB CM-250	Phenion® Fibroblast Culture Medium (250 mL)	\$ 60
K CM-250	Phenion® Keratinocyte Culture Medium (250 mL)	\$ 96
ALI CM-250	Phenion® Air-Liquid Interface Culture Medium (250 mL)	\$ 84
ALI CM PR-250	Phenion® Air-Liquid Interface Culture Medium w/o phenol-red (250 mL)	\$ 120
ALI CM HC-250	Phenion® Air-Liquid Interface Culture Medium w/o hydrocortisone (250 mL)	\$ 120
ALI CM AB-250	Phenion® Air-Liquid Interface Culture Medium w/o antibiotics (250 mL)	\$ 120
OS-REp ALI CM-125	Phenion® OS-REp Air-Liquid Interface Culture Medium (125 mL)	\$ 84
OS-REp ALI CM-250	Phenion® OS-REp Air-Liquid Interface Culture Medium (250 mL)	\$ 120
CS-3050	<b>NEW epiCS®</b> Culture Medium (50 mL)	\$ 36
CS-3051	<b>NEW epiCS®</b> Culture Medium (75mL)	\$ 54
CS-3052	<b>NEW epiCS®</b> Culture Medium (100 mL)	\$ 72
CS-3053	<b>NEW epiCS®</b> Culture Medium (125 mL)	\$ 90
CS-3030	<b>NEW epiCS®</b> MTT Assay Medium (25 mL)	\$ 18
CS-3031	<b>NEW epiCS®</b> MTT Assay Medium (50 mL)	\$ 36
CS-3151	<b>NEW epiCS®-M</b> Culture Medium (125 mL)	\$ 96
CS-3152	<b>NEW epiCS®-M</b> Culture Medium (500 mL)	\$ 312
CS-3130	<b>NEW epiCS®-M</b> MTT Assay Medium (25 mL)	\$ 18
CS-3131	<b>NEW epiCS®-M</b> MTT Assay Medium (50 mL)	\$ 36
<b>Accessories</b>		
FCA-15	Phenion® Franz Cell Adapter The tailored solution for permeation experiments with Phenion® FT Skin Models (fits to 15 mm PermeGear Franz Cells)	\$ 180
CS-5010	<b>NEW Nylon Meshes</b> (12 pieces)	\$ 84
CS-5011	<b>NEW Filter discs</b> (12 pieces)	\$ 84

We allow discounts in case of exceeding following order volumes:  
> \$ 3.600 = 5%, > \$ 7.200 = 7%, > \$ 12.000 = 9%

Contact for further information & orders:  
[www.phenion.com](http://www.phenion.com)

✉ [info@phenion.com](mailto:info@phenion.com)

